FR

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EN

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This Instructions for Use (IFU) is related to ALCIS wrench from range 2023-CL. It is an additional part of the IFU associated to ALCIS depth electrodes from the range 2069-EPC and associated accessories (ranges 2023 and 2036).

1. Products description

Performances:

- **Wrench REF 2023-CL and 2023-CL-C**
 - Ability to be assembled with guidance screws to allow them to be unlocked on the patient's skull.

Technical characteristics of the wrench:

Range 2023-CL  Year of CE marking: 2002	Wrench 	 	
		REF 2023-CL-C <i>Technical form 2023-AVG-C-FT</i>	Useful opening 3.5mm; to be used exclusively with ALCIS short guidance screws range 2023-VG-C .
		REF 2023-CL <i>Technical form 2023-APR-FT</i>	Useful opening 1.9mm to be used exclusively with standard guidance screws range 2023-VG .

2. Intended use

Wrench is intended to unlock the guidance screw before unscrewing by the screwdriver. As the wrench is an accessory of Depth Coagulation Electrode, (2069-EPC) It is used in the context of pre-surgical diagnosis of drug-resistant epilepsies.

3. Contraindications

No specific contraindications.

4. Recommendations / Warnings

The use of 2023-CL Wrench must be carried out by an experienced neurosurgeon and according to a controlled stereotactic procedure. The choice of the model is responsibility of the neurosurgeon. The patient must be placed under a continuous surveillance and in an adapted environment to avoid every risk of incidents.

- Do not “force” for accessories assembly or disassembly (assemblies must be performed without excessive stress)
- Check the products integrity and condition before each use.

Storage and maintenance

- Store in a cold dry place and away from sunshine at room temperature.
- Inspect packaging before use to detect any damage.
- Do not use in case of wear and/or damage of any kind. (risk of fall which can cause a deterioration or a breakage of the device)

Cleaning, disinfection and sterilization

Reusable accessories

any devices identified with any of the symbols shown here

(Refer to product label).

- Reusable devices are provided unsterile
- Before use, clean and sterilize according to a controlled and validated process (alkaline detergent and sterilization in an autoclave at 134°C for 18 minutes – prion cycle).

any devices identified with symbol shown here (refer to table 1 or product label)

- During the cleaning, take care to biological residues absence.
- After use, clean the reusable ancillaries (alkaline detergent).

Other important user information

- Check if the device is sold with this Instructions for uses to use it in recommended conditions.

Disposal procedure

There are infectious risks related to wrench disposal, according to WHO waste classification. To minimize the risks, the most common treatments for infectious waste are: steam treatment (autoclaving), integrated steam treatment (hybrid autoclaves or advanced steam treatment), microwave treatment, dry heat treatment, chemical treatment (internal shredding or chemical disinfectants) and incineration (combustion, pyrolysis or gasification). One of this method shall be followed to reduce the infectious risk. One of these methods should be followed to reduce the risk of infection.

5. Main complications

Any serious incident that as occurred in relation to the device should be reported to the manufacturer and the competent authority of the member state in which the user and/or patient is established.

However, no specific complications or potential adverse effects are identified for the wrench.

Main complications of interventional technics in stereotactic conditions still apply.

6. Warranty limit

ALCIS products are submitted to control tests at several stages of manufacturing process to ensure compliance with specifications. ALCIS can never be held responsible for any expenses or damages directly or indirectly resulting from:

- The non-respect of instructions, recommendations and warnings of the instruction for use,
- The use of accessories or devices unsuited or defective

7. Use procedure

Electrodes removal

Warning: To do without delay at the end of the thermo-coagulation procedure

Use with guidance screws 2023-VG:

- Check the absence of abnormalities with a suitable technique.
- Unscrew the cap.
- Remove the electrode.
- Use a wrench (2023-CL) on flat parts to unlock the screw.
- Use a short screwdriver (2023-TO-020) to finish removal of the screw.
- Close the cutaneous hole after screw removal.
- Reiterate the above operations for each screw.
- Check the cerebral structures with a suitable technique.

Use with short guidance screws 2023-VG-C:

- Check the absence of abnormalities with a suitable technique.
- Remove the retaining ring and unscrew the cap using a suitable screwdriver (2023-TO-BC).
- Then pull the electrode to remove it.
- Unlock the screw with a wrench (**2023-CL-C**), then unscrew totally with a short screwdriver (2023-TO-C-020).
- Close the cutaneous hole.
- Reiterate the same operations for each electrode.
- Check the cerebral structures with a suitable technique.

References table

Reference	Material	Packaging
2023-CL	Stainless steel	Unitary non-sterile
2023-CL-C		

8. Definition of labelling and packaging symbols

Table 1 – Symbols definition				
 REF	Manufacturer's Catalog Reference	 Keep away from sunlight	 LOT	Batch No
	Store in a dry place	 Do not use if damaged		Sterilize to re-use validated for x <i>nn</i> cycles (here 50)
	Manufacturing date	 Period of use After Opening (for this device 24 months)		Read Instruction For Use
	Manufacturer	 Warning		Steam sterilize before use: 134°C – 18min
	International phone number	 MD		Medical Device

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